ventura flow

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EN INSTRUCTIONS FOR USE

1. DEFINITION

Ventura flow is a flowable, light-curing, radio-opaque microhybrid composite with extremely low shrinkage stress. The composite material meets the EN ISO 4049 requirements.

2. COMPOSITION

Inorganic filler (barium glass 53 % wt, 28 % by volume, particle size between 50 – 3000 nm), BisGMA, BisEMA, catalysts, inhibitors and pigments.

3. INDICATIONS

1. Restorations with minimally invasive preparation technique.
2. Restorations of small cavities and extended fissure sealing.
3. Alternate restorations for undercut cavities.
4. Restorations of class III – V including wedged shaped defects and cervical caries.
5. Repair of fillings, veneers and methacrylate-based temporary restorations.
6. First layer of fillings for Class I and II.
7. Interlocking of loosened teeth.
8. Adhesive attachment of indirect composite and ceramic restorations.
9. Adhesive attachment of lingual retainers.

4. CONTRAINDICATIONS

Opened pulp, pulpitis, known allergy to methacrylate.

5. SIDE EFFECTS

In individual cases, contact allergies have been described in the case of products of a similar composition. To avoid pulp reactions, it is advisable to cover the exposed dentine in the pulp area with suitable underfilling materials (preferably with calcium hydroxide). Note: **Ventura Flow** does not contain TEGDMA or HEMA.

6.INCOMPATIBILITY WITH OTHER MATERIALS

Avoid underfilling materials which may hinder polymerisation owing to their ingredients. All phenolic compounds, such as zinc oxide eugenol, fall under this category.

7.PROCESSING STAGES

7.1. Tooth cleaning

Brush the tooth to be treated and the teeth neighbouring it with a fluoride-free toothpaste. Clean interdental spaces using strips and dental floss if necessary.

7.2. Shade selection

Under natural light, determine the shade using the VITA® SHADE GUIDE.

7.3. Drying

Dry sufficiently, apply rubber dam.

7.4. Cavity preparation

Prepare the cavity in the usual way. Undercuts and bevelled margins are recommended to improve adhesion and the margin fit. Clean and dry the cavity.

7.5. Underfilling

In the case of preparations near pulp, it is advisable to protect the pulp using an underfilling (e.g. calcium hydroxide). Do not use zinc oxide eugenol!

7.6. Etching

Apply an etching gel (e.g. **Ventura etching gel**) to the enamel and dentine according to the manufacturer instructions.

7.7 Adhesion

Apply a light cure adhesive (e.g. **Ventura unibond 2**) according to the manufacturer instructions.

8. RESTORATION WITH VENTURA FLOW

Remove syringe cap and screw on luer-lock cannula. Apply material with constant pressure on the plunger. Do not retract the plunger after application, retraction may cause air to be sucked into the material, thereby deteriorating the material properties.

Gradually adapt **Ventura Flow** up to a maximum layer thickness of 2 mm.

Light-cure every layer for 40 seconds. When using metal matrices, after removing them expose filling to light again from lingual and buccal.

Warning: After each usage, the application tip must be removed and the syringe tightly sealed with the original cap.

The specified exposure times refer to halogen or LED light-curing devices with a minimum light intensity of 500 mW/cm² and a wavelength of 400 – 500 nm. They apply to a maximum layer thickness of 2 mm and to all shades. The necessary exposure time may vary depending on the light source and its instructions for use. If in doubt, check the lamp's light output and the necessary exposure time before operation in vitro. During polymerisation, an inhibition layer forms on the surface; this must not be touched or removed, if other composite layers are going to be applied.

Light power	≥500 mW/cm ²	≥800 mW/cm ²
Color opaque	40 sec	30 sec
All other colours	40 sec	20 sec

9. FINISHING AND POLISHING

Trim filling with 40 µm and 12 µm diamonds, points or burs. Polish to high gloss with discs or silicone points. Interproximal finishing is accomplished by fine grit finishing strips. The filling can be prepared and polished immediately after polymerisation.

10. STORAGE

Protect light-curing products from strong sources of light and heat. **Ventura Flow** was developed for use at room temperature (20°C - 25°C / 68°F - 77°F). Store at temperatures between 4°C - 28°C / 39°F - 82°F. If room temperature exceeds 28°C / 82°F storing in the refrigerator is recommended. Do not freeze!

11. BATCH NUMBER AND EXPIRY DATE

The batch number should be specified to identify products in case of enquiries. Products should no longer be used once the expiry date has elapsed.

12. PRECAUTIONARY MEASURES

Close containers after each use using the right lid. Keep out of reach of children. When using tips, ensure before use that the cap has been removed. Commercially available medical gloves do not provide protection against the sensitisation effect of methacrylates. If the product comes into contact with the glove, remove the glove and dispose of it, wash your hands with water and soap immediately and put on a new glove. In case of an allergic reaction, seek medical advice. When polishing or removing composites, it is recommended to always use a water cooling system and a good extraction system, to ventilate the dental laboratory frequently and to wear masks with high particle filtration efficiency for small particle sizes.

13. EMERGENCY PROCEDURES

Contact with the skin: wash with soap and water. Contact with the oral mucous membrane: rinse with large quantities of water. Contact with the eyes: wash immediately with large quantities of water and consult an ophthalmologist.

14. HYGIENE

Use the application instruments for only one patient. Clean and disinfect all other instruments after each use. To avoid contamination, dispense the materials at a distance from the patient.

15. WARRANTY

Our liability is restricted to the quality of our products. In the case of a product being of defective quality, only its value is replaced. We accept no liability for other damage, namely that caused by non-compliance with the instructions for use or other improper handling or unintended use of a product. Before using the products, the user should check whether they are suitable for the intended purpose. He alone assumes all of the risks associated with using the product and bears sole responsibility for any damage caused as a result of this.

16. ADDITIONAL NOTES

If the user and/or patient become aware of serious incidents connected with application the product, they are to be reported to the manufacturer and the responsible authorities of the country in which the user and/or patient resides.

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