ventura crystalsil

REFERENCE	DESCRIPTION
3275	VENTURA CRYSTALSIL (2x50 ml) + 12 MIXING TIPS

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DESCRIPCIÓN

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EN INSTRUCTIONS FOR USE:

DESCRIPTION

Ventura crystalsil is a transparent material based on vinylpolisiloxanes suitable for all types of bite registrations that require transparent matrices.

INDICATIONS

Accurate occlusal registrations.
Key material for needle point registration.
Production of small model segments.
Silicone key for composite restorations.

DIRECTIONS OF USE

Insert cartridge onto dispensing gun following the manufacturer's instructions and remove cap. Extrude impression material until it exits both chambers at the same rate and then mount the appropriate mixing tip. An intra-oral tip may be placed on the mixing tip.

APPLICATION

Apply ventura crystalsil onto a transparent appropriate impression tray keeping the end of the mixing tip immersed in the material to avoid the formation of air bubbles and place in mouth. Carefully remove the set matrix after the recommended time. Rinse with water and blow dry with air. Prepare teeth according to the desired technique and following the manufacturer's instructions of use of the composite material. Apply composite onto the transparent matrix and position in mouth. Light cure through the matrix following the corresponding instructions of use. After use keep mixing tip on cartridge as it serves as a cap. Intraoral tips and mixing tips are for single use only.

DISINFECTING

Once bite registration is removed from mouth, rinse under cold water. It can then be disinfected with appropriate disinfecting solutions or sprays.

PRECAUTIONS

Contact with latex gloves may impair the setting reaction. Vinyl gloves are recommended.
Do not combine with condensation curing silicones.
Avoid contact of materials with eyes. In case of contact rinse eyes with plenty of water and consult a physician immediately.
Avoid spillage of materials on clothing as they cannot be removed chemically.
Do not use after the expiry date.

STORAGE

Store material at 10 - 25°C (50 - 77°F). Store bite registration at room temperature (max. 25°C/77°F).

LOT NUMBER AND EXPIRY DATE

The above are printed on each package and should be quoted in all correspondence which requires identification of the product.

TECHNICAL CHARACTERISTICS

PARAMETER	SPECIFICATION
Setting Time	Fast
Total working time*	1'00"
Time in mouth	2'00"
Total setting time*	3'00"
Dimensional change (%)	<0,3
Shore A hardness	70

* Attention: *Application is at 23°C ± 2 / 73°F ± 4, 50 ± 5% relative humidity. Higher temperatures reduce working and setting times and lower temperatures prolong them.

FOR PROFESSIONAL USE ONLY

These materials have been developed for use exclusively by dentists and must be handled according to the instructions of use. The recommendations included in the instructions correspond to the actual state of knowledge regarding dental techniques. The manufacturer cannot be held responsible for damages caused by other uses or incorrect handling. Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

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MADESPA, S.A. - Río Torviscal, 8 - 45007 TOLEDO - ESPAÑA (SPAIN)