AKRYON

- Ⓔ [EN](#) [Alginate dental impression material](#)
- Ⓔ [ES](#) [Material de impresión dental de alginato](#)
- Ⓔ [PT](#) [Material de impressão dentária de alginato](#)
- Ⓔ [IT](#) [Materiale per impronte dentali in alginato](#)
- Ⓔ [FR](#) [Matériau d'empreinte dentaire en alginate](#)
- Ⓔ [PL](#) [Materiał wyciskowy alginatowy do zastosowań stomatologicznych](#)
- Ⓔ [DE](#) [Alginate-Abformmaterial für Zahnmedizin](#)
- Ⓔ [NL](#) [Tandheelkundig afdrukmetaal op basis van alginaat](#)
- Ⓔ [LT](#) [Alginate odontologinė atspaudų medžiaga](#)
- Ⓔ [EL](#) [Οδοντιατρικό υλικό αποτύπωσης αλγινικού](#)
- Ⓔ [LV](#) [Algīnāta zobārstniecības nospiedumu materiāls](#)
- Ⓔ [BG](#) [Алгинатен стоматологичен отпечатъчен материал](#)
- Ⓔ [SK](#) [Alginátový dentálny odtlačkový materiál](#)
- Ⓔ [CS](#) [Alginátový dentální otiskovací materiál](#)
- Ⓔ [HU](#) [Alginát alapú fogászati lenyomatanyag](#)
- Ⓔ [FI](#) [Algiinipohjainen hammasjäljennösaine](#)
- Ⓔ [DA](#) [Alginatebaseret tandaftryksmateriale](#)
- Ⓔ [SL](#) [Alginate zobozdravstveni odtisni material](#)
- Ⓔ [SV](#) [Alginerat dentalt avtrycksmaterial](#)
- Ⓔ [HR](#) [Alginate stomatološki otisni materijal](#)
- Ⓔ [NO](#) [Alginatebasert dentalt avtrykksmateriale](#)
- Ⓔ [RO](#) [Material de amprentă dentară pe bază de alginat](#)
- Ⓔ [TR](#) [Aljinat bazlı diş ölçü materyali](#)

DESCRIPTION

AKRYON is a non-active medical device consisting of an elastic alginate-based dental impression material intended for obtaining negative impressions of oral structures in dentistry and dental prosthetics, through transient invasive use via the oral cavity.

The product is supplied as a powder which, when mixed with water according to the proportions indicated by the manufacturer, forms a mouldable paste capable of adapting to the morphology of the oral cavity and accurately reproducing anatomical details. After setting, the material forms an elastic and irreversible mould suitable for complete removal from the oral cavity.

INTENDED PURPOSE

AKRYON is an alginate-based dental impression material intended for obtaining negative impressions of oral structures.

Impressions obtained with AKRYON are used for the production of dental models, which are subsequently used for clinical evaluation, indirect diagnosis and the indirect fabrication of dental devices.

The product is designed for professional use in the dental and prosthetic field, and is used for the time required to take the impression and removed immediately after setting. AKRYON does not in itself provide a diagnostic result, but acts as an auxiliary tool within the clinical and technical process of planning and manufacturing dental devices.

CLINICAL BENEFITS

Impression material intended to accurately record the dimensions of oral tissues or prosthetic components and the spatial relationships between them.

INTENDED USER / FOR PROFESSIONAL USE ONLY

The product under evaluation is intended to be used by qualified healthcare professionals in the fields of dentistry and dental prosthetics.

Preparation of the material, impression taking and removal of the impression must be carried out by personnel with the appropriate knowledge and experience to ensure accurate reproduction of the oral structures and patient safety, and in accordance with the manufacturer's instructions for use.

The product is not intended for use by patients or for self-application.

The intended user must follow the manufacturer's instructions regarding the mixing ratio, working and setting times, as well as the conditions of use and removal of the material, in order to ensure an appropriate result and minimise the risk of impression errors.

Use in specific populations, such as paediatric or geriatric patients, is subject to the clinical judgement of the responsible healthcare professional.

WARNINGS AND PRECAUTIONS

The product must not be used under conditions other than those described in the manufacturer's instructions for use.

Handle the powder carefully. Do not inhale the powder and avoid contact with the eyes. The product must not be used in the neonatal population.

All products must be used at room temperature. Allow the product to reach room temperature before use. High temperatures reduce (accelerate) the working and setting times, whereas low temperatures increase (slow down) these times.

Do not ingest. Keep the container tightly closed.

Do not use if the powder is damp or has formed lumps. Do not reuse the impression.

CONTRAINDICATIONS

Do not use in patients with a history of sensitivity and/or allergic reactions to any of the components of the product.

Do not use in the presence of extensive lesions, active ulcerations, open wounds or significant bleeding in the oral cavity; consider an alternative where appropriate.

The product must not be used in the neonatal population.

DIRECTIONS FOR USE

1. Preparation and mixing

AKRYON is supplied in powder form and must be mixed with clean potable water in the powder/water ratio specified by the manufacturer.

The material must be mixed for the recommended mixing time until a homogeneous, lump-free paste with a consistency suitable for clinical use is obtained.

Preparation must be carried out in accordance with the environmental conditions specified in the Instructions for Use, in order to ensure correct material performance during the working and setting phases.

2. Impression taking

Once the mixture has been prepared, the alginate paste must be loaded into a suitable impression tray and applied in the patient's oral cavity within the specified working time. During this phase, the material adapts plastically to the morphology of the oral structures, allowing accurate reproduction of the anatomical details.

The alginate remains in temporary contact with the oral cavity until the setting process is complete.

3. Setting and removal

Setting of the material takes place through an irreversible physicochemical process, after which an elastic and dimensionally stable gel is obtained.

Once complete setting has been achieved, the impression must be removed from the oral cavity in one piece, avoiding any deformation that could compromise the quality of the anatomical record.

4. Disinfection

Once the impression has been removed from the oral cavity, it must be carefully rinsed for a few seconds under running water at room temperature (do not immerse). It must then be carefully dried before proceeding with the corresponding disinfection.

Disinfection may be carried out either by spraying with disinfectants specifically intended for this product or similar products, or by immersing it in a specific disinfectant solution. In both cases, the specific recommendations of each manufacturer must always be followed.


5. Use of the impression and subsequent handling


The impression obtained with AKRYON is intended for immediate or short-term use for pouring with dental plaster in order to obtain a positive model.

The material is compatible with dental plasters and does not inhibit their setting nor affect the accuracy of the model obtained, provided that the manufacturer's instructions are followed.

6. Disposal after use

Following impression taking and its intended use, the material must be regarded as potentially contaminated, as it has been in contact with the oral cavity and biological fluids.

CLINICAL USE SPECIFICATIONS

Ratio (powder:water)	1:2
Mixing time	30"
Working time	1'30"
Setting time	3'30"
Time in contact with the mucosa	1'30"
Total time	3'30"
Shelf life	3 años

STORAGE CONDITIONS

The recommended storage temperature is between 5 °C and 25 °C.

STORAGE CONDITIONS AFTER OPENING

Keep in a dry place and in an airtight container.

SYMBOLS

	Manufacturer
	Use-by date
	Batch code
	Do not use if the package is damaged and consult the instructions for use
	Keep dry
	Temperature limits for storage
	Consult the instructions for use
	Medical device
	Unique Device Identifier

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MADESPA, S.A.