ventura **Impress 90** *lab*

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INSTRUCTIONS FOR USE

1. DESCRIPTION

VENTURA IMPRESS 90 is an elastic addition-curing silicone (polyvinyl siloxane) for use in dental laboratories.

It is supplied as two high-viscosity putties (Base + Catalyst), which are easily mouldable and mixed manually.

VENTURA IMPRESS 90 LAB is suitable for all applications where an accurate reproduction of a surface is required, such as registrations or moulds. It also ensures high detail reproduction under all conditions of use and, together with the high degree of hardness achieved by the polymerised material and its high elastic recovery, makes it particularly suitable for applications requiring greater rigidity and for all cases where the registration needs to be maintained for a long period of time.

Surfaces obtained with **VENTURA IMPRESS 90** can be easily scanned using the most common CAD/CAM/CIM systems.

The applications of **VENTURA IMPRESS 90** include:

- Occlusal registration keys.
- Matrices for the fabrication of:
 - Partial dentures.
 - Temporary restorations.
 - Bars.
 - Anterior teeth (immediate dentures).
 - Gingival reproduction (indirect technique).
 - Ceramic pressing.
- Control keys for metal frameworks.
- Control keys for diagnostic wax-ups.
- Keys for processing light-curing resins.
- Isolation of teeth in the flask.
- Support for denture repairs.
- General duplication applications.

2. INSTRUCTIONS FOR USE

Before use, properly prepare the surfaces with which **VENTURA IMPRESS 90** will come into contact. The surfaces must be dry and free from oily residues in order to obtain an accurate and well-detailed reproduction.

Then, using the supplied measuring scoops (using one colour for the Base and another for the Catalyst), take equal amounts of Base and Catalyst as required (mixing ratio 1:1) and manually mix the two components for 30 seconds until a completely homogeneous, streak-free colour is obtained.

Next, shape the material according to the work to be carried out, taking care not to exceed the working time specified in the technical characteristics table.

The setting time is approximately 5–6 minutes. Therefore, do not remove the product until this time has elapsed before proceeding with the subsequent processing stages.

3. TECHNICAL CHARACTERISTICS

	Ventura Impress 90
Appearance after mixing	Homogeneous blue paste
Mixing time	30"
Working time	2'30 – 3'00"
Setting time	5'00 – 6'00"
Shore A Hardness	> 88
Deformation under compression	0,6 - 1 %
Linear dimensional change	- 0,04 - +0 %

4. WARNINGS/PRECAUTIONS

The stated times are measured from the start of mixing at 23°C (73°F). Larger amounts of catalyst and/or higher material and/or room temperatures shorten the stated times, whereas smaller amounts of catalyst and/or lower temperatures extend them.

Wear protective gloves and protective clothing. **Use vinyl/nitrile gloves. Avoid contact between the silicone and latex gloves or components made from latex. If latex gloves have previously been used, wash your hands thoroughly. Latex inhibits the setting of polyvinylsiloxanes.**

5. STORAGE

Do not store at temperatures above 25°C (77°F). Protect from direct sunlight and do not use after the expiry date. Close the containers tightly after use, taking particular care not to interchange the lids of the two components or the measuring scoops, as cross-contamination between the two putties may trigger partial polymerisation, which could damage the product and render it unusable.

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MADESPA S.A.