

A close-up photograph of a green and red iguana, likely a Spiny-tailed Iguana, perched on a branch. The iguana's head is the central focus, showing its textured scales in shades of green, yellow, and red. Its eye is prominent, and its mouth is slightly open. A large, white, stylized letter 'M' is overlaid on the left side of the image, partially obscuring the iguana's body. The background is a soft, out-of-focus blue and green.

Madespa

Impression
material

Impress 2.0

Heavy body addition silicone

Especially indicated for:

Double impression technique, using two different viscosities of silicone.

Features and benefits:

Type mouldable putty, very high viscosity, more comfortable for both insertion and removal from the mouth.

Great resistance to deformation and tearing.

Excellent elastic recovery.

Greater dimensional stability.

Manual mixing system with base and catalyst in putty form.

Formats:

Putty Soft Regular set: 250 ml base + 250 ml catalyst

Putty Soft Fast set: 250 ml base + 250 ml catalyst

Properties:

	Putty Soft Regular set	Putty Soft Fast set
Working time:	1'45"	1'15"
Total setting time:	3'	2'
Shore A hardness:	68	68
Colour:	Yellow	Blue



Impress 2.0

Light body addition silicone

Especially indicated for:

Double impression technique, using two different viscosities of silicone.

Features and benefits:

Flowable silicone, low viscosity for precise detail reproduction.

Hydrocompatible: captures detail in a moist oral environment.

High thixotropy.

Excellent elastic recovery.

Great resistance to deformation and tearing.

Dual cartridge format with auto-mix tip: ensures a homogeneous mix, precise dosing and material saving.

Formats:

Light Body Regular set:

- 2x50 ml + 12 mixing tips
- 4x50 ml + 4 mixing tips

Light Body Fast set:

- 2x50 ml + 12 mixing tips
- 4x50 ml + 4 mixing tips

Properties:

	Light Body Regular set	Light Body Fast set
Working time:	2'	1'
Total setting time:	3'30"	2'
Shore A hardness:	56	56
Colour:	Green	Orange



Bite 2.0

Addition silicone for bite registration

Epecially indicated for:

Reliable bite registrations and without deformations.

Features and benefits:

It is a vinylpolysiloxane suitable for all types of occlusion techniques.

High thixotropy.

Easy to polish after setting.

Fast setting for patient comfort, only 45" in the mouth.

Auto-mix system: easy, comfortable, and bubble-free.

Formats:

Ventura Bite 2.0 (2 x 50 ml) + 12 mixing tips.

Properties:

Working time:	12"
Total setting time:	45"
Shore A hardness:	89
Colour:	Lavender



Addition silicone for bite registration for all types of occlusion techniques, easy to use and accurately reproduces occlusal details.

Crystalsil

Transparent silicone for matrices

Especially indicated for:

The creation of transparent matrixes for aesthetic restorations in light-curing composites, using direct and indirect techniques.

Features and benefits:

It faithfully reproduces the wax-up for light-curing composite restoration thanks to its perfect flow during application and high degree of transparency. It has high rigidity, preventing deformation.

Formats:

Ventura Crystalsil (2 x 50 ml) + 12 mixing tips.

Properties:

Working time:	1'
Intraoral time:	2'
Total setting time:	3'
Shore A hardness:	70 ± 2
Colour:	Transparent



Its great transparency and processing time make it particularly effective for discolouration, fractures, rotation, conical teeth and other types of dental malformations.

Top 2.0

Condensation silicone

Especially indicated for:

The double impression technique or simultaneous impression technique with two viscosities in fixed prostheses due to its high density after polymerization.

Features and benefits:

Condensation silicone offering high precision and an ideal balance between working time and setting time.

Top 2.0 Putty Yellow putty with a mint fragrance. Excellent hardness. Fast working and setting time.

Top 2.0 Light Blue fluid silicone. It provides a perfect color contrast when used with the putty. Great flow and precision of the corrective materials.

Formats:

Ventura Top 2.0 Putty - 900 ml
Ventura Top 2.0 Catalyst gel - 60 ml
Ventura Top 2.0 Light - 140 ml
Ventura Top 2.0 Kit

Properties:

	Top 2.0 Putty	Top 2.0 Light
Working time:	1'15"	1'30"
Intraoral time:	2'30"	4'
Total setting time:	3'45"	5'30"
Shore A hardness:	65 ± 2	45 ± 2
Colour:	Yellow	Blue



Akryon

Alginate for dental impression

Especially indicated for:

Dental impression procedures requiring ease of use, good detail reproduction and reliable results, such as study models, removable prostheses, antagonists, orthodontics and splints

Features and benefits:

Dental alginate designed to provide a homogeneous mix, easy handling and accurate detail reproduction.

It offers an optimal consistency to facilitate impression taking and to deliver reliable results in routine dental procedures.

Formats:

Akryon Alginate - 500 g bag

Properties:

Mixing ratio (powder:water):	1:2
Flavour:	Melon
Colour:	Green
Mixing time:	30"
Working time:	1'30"
Total setting time	3'30"





*M*adespa

Río Torviscal, 8 - 45007 - Toledo

España - Spain

Tel. (+34) 925.24.10.25

Fax. (+34) 925.24.10.28

E-mail: info@maadespa.com

international@maadespa.com

www.madespa.com

Rev.02 2026-09

